

**Equipment  
-25**

**Optical Microscope**

**Picture**



Leica MS5



Olympus BX51M

**Brand**

Leica  
Olympus

**Model**

MS5  
[BX51M](#)

**Specification**

The Leica MS5 stereo microscope offers a 5 steps magnification changer. This provides .63x, 1x, 1.6x, 2.5x, and 4x magnification stops. This microscope is ideal for situations where you need to know you are for certain at a specific magnification and not in-between magnifications like you could be with a stereozoom microscope. We offer this microscope with 10x eyepieces

Click link for Olympus BX51M model as above.

**Status**

Available

**Charge Rate (RM)**

| Unit | Category   | Per Hour |
|------|--|----------|
| A    | UPM's Student, Staff, Association, UPM's Faculties, Institution and Schools        | 105.00   |
| B    | Government Agencies, Statutory Body  | 135.00   |
| C    | Private sectors, Individual, UPM's Subsidiaries, Union, Government's Business Body | 170.00   |

## Contact

|                         |   |
|-------------------------|---|
| <b>Person In charge</b> | Muhammad Wildan Ilyas Mohamed Ghazali   |
| <b>Tel. No</b>          | 03-9769 4681  |
| <b>E-mail</b>           | wildan@upm.edu.my   |
| <b>Makmal</b>           | Makmal Kekuatan Bahan, KMP 008  |
| <b>Location</b>         | <a href="https://goo.gl/maps/P6ukENrKnDEF32rJ8">https://goo.gl/maps/P6ukENrKnDEF32rJ8</a> |